





Cassette Flat panel detector

Wired 43×43CM



4343ZCV

- Stable
- High-efficiency
- Innovative



Specification parameter	
TECHNOLOGY	a-Si
Scintillator	CsI
Effective Imaging Area (mm)	43x43cm
Active Matrix	3072×3072
Effective Pixel Pitch (μm)	139
A/D Conversion (bits)	16
Spatial Resolution	3.6LP/mm
X-ray working Range	40-150KV
Image Acquisition Time (WIRED)	1s
Data interface	GigE
Adapter output	DC 24V,60W
Adapter input	AC100-240V,50-60Hz
Power	20W
Dimensions	46.0x46.0x1.5cm
Weight	3.5kg
Housing material	Carbon, Plastic
Waterproofing grade	IP54
Working environment	5-35°C,10-75% RH
Schematic Drawing	

